



## Oxford AstraZeneca ChAdOx1 nCoV-19

### From the Green Book:

*"The AstraZeneca vaccine should be stored at +2°C to +8°C and has a shelf life of 6 months. The vaccine does not contain any preservative. After first opening the vial, it should be used within 6 hours. The vaccine may be stored between 2°C and 25°C during this period. After this time, the vial must be discarded."*

### Easy, Low Cost Monitoring: EL-USB-VAC

The EL-USB-VAC data logger is perfect for monitoring individual fridges set to between +2°C and +8°C.

- Highly accurate glycol bottle thermistor probe that mimics the temperature response of individual vaccine containers.
- Coloured indicators show the logger is working correctly, and whether an alarm has been triggered.
- Battery included giving 2 years of uninterrupted service.
- Calibration certificate supplied to ensure accuracy.

*To check your measured data in detail, just plug the logger directly into your PC or laptop, and use our free EasyLog software to graph, analyse and download the readings.*

Available NOW for just

**£135.95 Each**



### WiFi Logger: EL-WiFi-VAC

Get critical vaccine condition monitoring text alerts straight to your pocket with our vaccine logger SMS functionality.

Cloud-connected wireless data loggers take just minutes to set up, providing fully automated monitoring across different storage areas or even different buildings. Your data is available 24/7 from any internet-connected device, and instant notifications via email and SMS messages mean you and your team are aware the moment a problem occurs.

- Use your existing wireless network, data is encrypted and secure.
- Rechargeable internal battery pack, and can also be permanently powered.
- Supplied with glycol probe and calibration certificate.

Available NOW for just

**£169.95 Each**

*To simultaneously monitor two fridges or storage zones, the EL-WiFi-VAC2 connects two probes from one data logger.*

 The EasyLog Cloud is FREE for one device, and costs as little as £0.38 per month for additional loggers.



To order your loggers call us on **01794 884 567** or visit **lascarelectronics.com** and use the search bar at the top of the page to look for EL-USB-VAC or EL-WiFi-VAC. Please also call us for technical help or to learn more about our other vaccine monitoring products.

Please turn over for information for the Pfizer and Moderna vaccines.



## Pfizer BioNTech BNT162b2 (Also suitable for Moderna mRNA-1273)

### From the Green Book:

*"The Pfizer BioNTech COVID-19 mRNA Vaccine BNT162b2 should be stored in a freezer at -80°C to -60°C (-90°C to -60°C in the thermal container). Shelf life is 6 months at -80°C to -60°C. Frozen vials should be transferred to 2°C to 8°C to thaw; a 195 vial pack may take 3 hours to thaw. Alternatively, frozen vials may also be thawed for 30 minutes at temperatures up to 25°C for immediate use. After thawing, stability data have demonstrated that undiluted vaccine can be stored for up to 5 days at 2°C to 8°C. Once thawed, the vaccine cannot be re-frozen."*

## Simple, Low Cost Monitoring: EL-USB-ULT-LCD

Available NOW for just

**£84.95 Each**

The EL-USB-ULT-LCD data logger is designed for monitoring vaccines in dry ice storage, down to -100°C.

- Highly accurate temperature probe can reach inside insulated containers.
- Coloured indicators show the logger is working correctly, and whether an alarm has been triggered.
- Battery included giving 2 years of uninterrupted service.

*To check your measured data in detail, just plug the logger directly into your PC or laptop, and use our free EasyLog software to graph, analyse and download the readings.*



## WiFi Logger: EL-WiFi-ULT

Get critical vaccine condition monitoring text alerts straight to your pocket with our vaccine logger SMS functionality.

Cloud-connected wireless data loggers take just minutes to set up, providing fully automated monitoring across different storage areas or even different buildings. Your data is available 24/7 from any internet-connected device, and instant notifications via email and SMS messages mean you and your team are aware the moment a problem occurs.

- Measures down to -100°C.
- Use your existing wireless network, data is encrypted and secure.
- Rechargeable internal battery pack, and can also be permanently powered.
- Supplied with glycol probe and calibration certificate.

*To simultaneously monitor two fridges or storage zones, the EL-WiFi-DULT connects two probes from one data logger.*



The EasyLog Cloud is FREE for one device, and costs as little as £0.38 per month for additional loggers.



Available NOW for just

**£139.95 Each**

To order your loggers call us on **01794 884 567** or visit [lascarelectronics.com](http://lascarelectronics.com) and use the search bar at the top of the page to look for EL-USB-VAC or EL-WiFi-VAC. Please also call us for technical help or to learn more about our other vaccine monitoring products.

Please turn over for information for the Oxford AstraZeneca vaccine.